



# BIOLOGICAL INDICATOR I N C U B A T O R



## Biological Indicator Incubator LBII-100 Series

Labocon Biological Indicator Incubator LBII-100 Series is designed for evaluating the sterility of instruments, vials, media, liquid cycles and other products in short span of time. It provides easy verification of the incubation temperature.

### Features

- Designed for incubation of EO biological indicators
- Customized wells for placement of biological indicator vials
- Easy to set incubation temperature
- Aluminium blocks available to meet experimental requirement

### Application

Labocon Biological Indicator Incubator is used for sterilization processes in food and beverage industries, medical and healthcare sectors, and pharmaceutical manufacturing applications.



LBII-101



LBII-102

## Specification

Model	LBII-101	LBII-102
Temp. Range	RT +5-80°C	RT +5-100°C
Temp. Accuracy	0.35°C at 57°C	
Temp. Uniformity	0.5°C	
Temp. Resolution	0.1°C	
Time for Heating	≤6min (20°C to 57°C)	
Timer Range	1min-99h59min	
Display	LCD	
Overall Dimension	110x150x80 mm	120x152x112 mm
Power	120W	
Power Supply	240V, 60Hz	
Weight (Net/Gross)	0.5/1 kg	0.85/1.4 kg
Catalog No.	9341239105	9341240105

## Accessories

Model	Description
LB-83	15xφ10mm
LB-84	15xφ12mm
LB-85	6xφ15mm
LB-86	6xφ16mm
LB-87	24xφ8.8mm
LB-88	18xφ9mm



**LABOCON SCIENTIFIC LIMITED**

18A Melton Road, Leicester, LE4 5EA, UNITED KINGDOM

 +44 203 3724877 |  [info@labocon.com](mailto:info@labocon.com) |  [www.labocon.com](http://www.labocon.com)