



— — — — — **MINI** — — — — —  
**VACUUM PUMP**

### MINI VACUUM PUMP LMVP-100 SERIES

Labocon Mini Vacuum Pump LMVP-100 Series uses big starting torque design and is especially designed for low temperature working environment. It has a compact design and requires less storage place.

#### Features

- ▶ Compact, reliable and powerful
- ▶ Especially designed for low temperature working environment
- ▶ Space-saving and ideal for systems integration
- ▶ High pumping speeds with small dimensions
- ▶ Eco-friendly exhaust gas free of oil mists
- ▶ Quiet, low-vibration operation
- ▶ Easy to maintain

#### Application

Labocon Mini Vacuum Pump can be used for vacuum packaging, vacuum freezing, gas analysis and other vacuum evacuations.



LMVP-101

## Specification

Model	LMVP-101	
Pumping Speed	28.31 L/min	22.65 L/min
Rotating Speed	1440 rpm	1720 rpm
Ultimate Pressure	5 Pa	
Oil Capacity	220 ml	
Dimension	243x114x207 mm	
Power	90 W	120 W
Power Supply	110 V, 60 Hz	220 V, 50 Hz
Weight	6.5 kg	
Catalog No.	9204122542	



### LABOCON SYSTEMS LIMITED

Fowler Avenue, The Hub, Farnborough Business Park  
Farnborough, GU14 7JF, United Kingdom

+44 203 3724850 | 
 info@labocon.com | 
 www.labocon.com