




PLATELET SHAKER

PLATELET SHAKER LPS-100 SERIES

Labocon Platelet Shaker LPS-100 Series is especially designed to meet the demands of various clinical laboratories. They are made of high-grade steel with stainless steel trays for easy cleaning and durability. It is supplied with multiple trays to accommodate 12-48 platelet bags and is low maintenance with quiet operation.

Features

- Platelet bags can be loaded easily on the stainless steel drawer
- Mesh design provides air circulation
- Exhibits a steady and smooth 360° movement
- Stainless steel shelves can be easily cleaned
- Easy monitoring of speed
- Equipped with a tachometer
- Optional: Basket



Application

Labocon Platelet Shaker is widely used in blood banks and clinical research laboratories.

Specification

Model	LPS-101	LPS-102	LPS-103	LPS-104
Capacity	12 bags	24 bags	36 bags	48 bags
Shaking Speed Range	0-9 rpm			
Shelves	1	2	3	4
Overall Dimension	370x255x290 mm	510x255x290 mm	650x255x290 mm	790x255x290 mm
Power	30 W			
Power Supply	220 V, 60 Hz			
Weight	7.5/8.5 kg	9/10 kg	10/11 kg	11/13 kg
Catalog No.	9173107162	9173108162	9173109162	9173110162

PLATELET SHAKER LPS-200 SERIES

Labocon Platelet Shaker LPS-200 Series is especially designed to meet the demands of various clinical laboratories. They are made of high-grade steel with stainless steel trays for easy cleaning and durability. It is supplied with multiple trays to accommodate 42-84 platelet bags and is low maintenance with quiet operation.

Features

- Platelet bags can be loaded easily on the stainless steel basket
- Mesh design provides air circulation
- Side-by-side 38mm movement of the storage frame
- Exhibits a steady and smooth 360° movement
- Equipped with a ball bearing drive motor and an internal fan



Application

Labocon Platelet shaker is widely used in blood banks and clinical laboratories.

Specification

Model	LPS-201	LPS-202
Capacity	42 bags	84 bags
Shaking Speed Range	70 rpm	
Shelves	7	
Overall Dimension	460x350x360 mm	840x350x360 mm
Power	30 W	50 W
Power Supply	220 V, 60 Hz	
Weight	26/28 kg	39/49 kg
Catalog No.	9173111162	9173112162



LABOCON SYSTEMS LIMITED

Fowler Avenue, The Hub, Farnborough Business Park
Farnborough, GU14 7JF, United Kingdom

 +44 203 3724850 |  info@labocon.com |  www.labocon.com