



Three Hole Water Bath

LABOCON SYSTEMS

3rd Floor, 207 Regent Street, London, W1B 3HH, UK

Email: info@labocon.com, Web: www.labocon.com

Three Hole Water Bath LWTB-100 Series



Labocon LWTB-100 Series of Three Hole Water bath is an economical water bath with three types of temperature settings for low, medium and high temperature. It is equipped with a drainage valve, thereby making it more convenient to use.

Features:

- * Digital PID controller
- * Digital display
- * Internal material made up of seamless stainless steel
- * External material made up of powder coated steel plate
- * Stainless steel lid



Application:

Labocon Three Hole Water Bath LWTB-100 Series is extensively used for warming reagents, serological research, protein analysis, molecular biology, genetic studies

Specification:

Model	LWTB-101	LWTB-102	LWTB-103
Capacity	6 L x 3	11 L x 3	22 L x 3
No. of Holes	3		
Temperature Range	Ambient +5°C to 50°C (Low Temperature) 40°C to 70°C (Medium Temperature) 60°C to 99°C (High Temperature)		
Timer Range	0-99 hr 59 min		
Internal Dimension	(155x300x150) x 3 mm	(240x300x150) x 3 mm	(259x500x50) x 3 mm
Overall Dimension	690x470x360 mm	960x470x360 mm	1050x670x360 mm
Power	1400 W	2300 W	4500 W
Catalog No.	9501371470	9501372470	9501373470